
METHODE BRIGITTE KETTNER

Thank you for purchasing the Sol Eye Massager. Please read this instruction manual carefully before using this product and keep it for future reference. This tool promotes collagen regeneration, blood circulation, helps relieve eye fatigue. It improves dark circles and under eye bags. Can be used on lips lines and nasolabial folds.

CHARGING YOUR DEVICE:

Plug your device into a USB port with the provided USB cable and charge for 2 hours. During the charging process, the display screen battery symbol flashes. The battery symbol stops flashing after charging is completed. (Make sure to remove the protective film from the display screen area)

DEVICE FEATURES:

Metal massage head with red LED light
LCD Screen
USB charging port
Upshift button
Downshift button
On/Off key
Mode selection button

Red				Light: 630nm-635nm
High	Frequency	Vibration:	10,000/per	minute
Adjustable Heat: 37° C -45° C/ 98.6° F- 113° F				

WORKING MODES

Mode 1: Vibration + Heat + Red Light

Using 10,000 times / minute high-frequency vibration massage, 37 ° C -45 ° C heat and red LED to enhance cell activity, promote cell regeneration and repair. Triple-effect, brightening and skin rejuvenation, helping decrease wrinkles, increases the elasticity of the eye muscles.

Mode 2: Heat + Red Light

The heat (37 ° C -45 ° C) wakes up sleeping eye muscles, and the LED red light enhances cell activity, promotes collagen regeneration, and repairs damaged cells.

Mode 3: Vibration + Heat

10,000 times / minute high-frequency vibration massage and hot compress (37 ° C -45 ° C) stimulates eye muscles and is great for product penetration.

USAGE:

1. Take an appropriate amount of eye cream and spread it evenly around the eye.
2. Hold the on/off button down for 3 seconds to turn on device and short press (on/off button) to switch the mode, and press "+" / "-" to adjust the temperature as needed.
3. Hold the tool so the metal tip is facing you or facing the client if being done in treatment.

4. Use upward and outward motions around the eye. For lower eye puffiness start on the inner eye and move out and downward for lymphatic drainage.
4. To turn off hold the on/off button down for 3 seconds. Wipe the massage head with an anti-bacterial wipe.

MAINTENANCE

1. Please apply the skin care products evenly before each treatment to avoid excessive skin care products accumulating in the massage head.
2. Please wipe the product clean with a dry cloth after each treatment is completed.
3. Please don't pour water on the device or submerge device in water.
4. Please don't use corrosive to clean the instrument.
5. Please store it in a cool dry place.
6. Please don't use the device if it is damaged.

CAUTIONS:

1. This product should not be used in the following: children, pregnant women, patients with epilepsy, tumors and cancer patients, patients with cardiac pacemakers, implanted with cardiac defibrillators, and other devices with implanted in vivo. If you have any medical questions, such as serious illness, epilepsy, etc., please consult your doctor before using.
2. If your eyes have inflammation, such as conjunctivitis, do not use this product. If you feel abnormalities in your eyes, eyelids, or skin during use or after use, please consult your medical specialist.
3. Remove all metals, such as earrings, necklaces, bracelets etc. before use.
4. Do not use while charging.
5. If using on sensitive skin, please test an area first.
6. Do not repair, disassemble, or modify this product.

WARRANTY:

All MBK products have a 2-year warranty. If you have any issues with your device, please contact MBK at laura.porter@telefieldna.com.